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CellColt Care

REF

PSH1DXXYY.C
PSH1RDXXYY.C
PSH1DXXYY.C
PSH1RD2XXYY.C

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ENGLISH

Indications and Description

1. Before use, check that the package is intact and not damaged.
2. Check on the label that the length and diameter of the needle correspond to what you want.
3. Open the package and remove the needle protection.
4. Disinfect the area intended for introduction and possibly perform local anesthesia. To facilitate the insertion of the needle, it is possible to make an incision on the skin using a scalpel.
5. With the trocar stylet inserted, introduce the needle through the skin until it reaches the cortical bone (periosteum).
6. Applying light but firm pressure, advance the tip of the needle through the periosteum inside the cortex, rotating the needle in an alternating clockwise and counterclockwise motion.
7. Once the cortical membrane has been passed (a significant decrease in resistance to penetration into the tissue is anticipated):
 - **Versión without blunt stylet (PSH1Dxxyy)**: go to step 8
 - **Versión with blunt stylet (PSH1D2xxyy)** extract the trocar stylet by turning the stylet cap counterclockwise. [Fig2] Insert the blunt mandrel into the hole in the handle by hooking it to the handle with a clockwise rotation.
8. Position the stopper 3 centimeters from the skin (corresponding to 3 centimeter marks). The stopper has the function of predetermining the correct insertion depth of the cannula so that all the holes go beyond the cortex. [Fig3]
9. Advance the needle with alternating clockwise and counterclockwise movements and continue until the stopper reaches the skin.
10. Remove the blunt mandrel by turning the mandrel cap counterclockwise.
11. Open the sliding block laterally so that it protrudes from the side of the handle. [Fig4]
12. Insert the internal cannula with an orange luer attachment into the hole in the handle and press it until it stops, making sure that the arrow positioned on the orange handle is inserted into the handle on the side of the MDL logo.
13. Connect a luer lock syringe to the Luer connection on the handle of the internal cannula.
14. Aspirate by pulling back the plunger of the syringe. **ATTENTION:** keep the syringe steady during aspiration to avoid inadvertently extracting the cannula from the handle. Suction about half of the fluid you intend to withdraw.
15. Stop suction, press the sliding block towards the handle [Fig5-Fig6] in order to cause the internal cannula to retract. In this way the suction holes used in the first phase are closed and others are opened in different positions. Make sure that in the new position the cannula always remains at the limit switch and in contact with the sliding block.
16. Repeat the aspiration by pulling back the plunger of the syringe until the desired quantity is reached.
17. Unscrew the syringe and remove it with the aspirated sample.
18. Insert the luer cap on the luer connection placed on the orange handle of the internal cannula in order to avoid blood leakage during needle extraction.
19. Extract the needle from the bone with light rotational movements.
20. After aspiration, it is indicated for use in gynecological, tissue, plastic, reconstructive, regenerative, reparative and scarring procedures and treatments
21. If the needle is equipped with NeedleVISE®, place the NeedleVISE® on a stable surface and insert the tip of the needle into the NeedleVISE®.
22. Discard the device according to specific hospital procedures.

Warnings and Precautions

For use by medical specialists only. The use procedure provided is for information purposes, the professional must evaluate the methods of use based on his own clinical experience and the type of intervention to be carried out; any changes are at the discretion of the specialist. Check the expiration date before use. The medical devices used are at high risk of contamination. Pay close attention to avoid accidental punctures. Treat the medical device as contaminated hospital waste. Dispose of the product in compliance with the laws and regulations in force. Do not use for an application other than that of soft tissue biopsy.

DO NOT USE THE DEVICE UNDER MRI (Magnetic Resonance Imaging).

This instrument has been manufactured for single use only. Reuse can cause loss of mechanical integrity and affect the functionality of the device. The M.D.L. SRL declines any responsibility for the reuse of the device.

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